

Product Datasheet

Mouse CD31 ELISA Kit (orb409220)

Description	Mouse soluble platelet endothelial cell adhesion molecule 1,sPECAM-1 ELISA Kit
Reactivity	Mouse
Range	3.12 pg/mL-200 pg/mL
Target	PECAM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma
Uniprot ID	Q08481
Sensitivity	0.78 pg/mL
Expiration Date	Please enquire.