

Product Datasheet

Mouse CSF3 ELISA Kit (orb409231)

Description	Mouse Granulocyte Colony Stimulating Factor,G-CSF ELISA Kit
Reactivity	Mouse
Range	7.8 pg/mL-500 pg/mL
Target	CSF3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma,tissue homogenates
Uniprot ID	P09920
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.