

Product Datasheet

Mouse CSF2 ELISA Kit (orb409232)

Description	Mouse Granulocyte-Macrophage Colony Stimulating Factor,GM-CSF ELISA Kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	CSF2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P01587
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.