

## Product Datasheet

### Mouse XCL1 ELISA Kit (orb409237)

|                          |                                          |
|--------------------------|------------------------------------------|
| <b>Description</b>       | Mouse lymphotactin,Lptn/LTN<br>ELISA Kit |
| <b>Reactivity</b>        | Mouse                                    |
| <b>Range</b>             | 0.156 ng/ml - 10 ng/ml                   |
| <b>Target</b>            | XCL1                                     |
| <b>Note</b>              | For research use only                    |
| <b>Application notes</b> | Detection Wavelength: 450 nm             |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates     |
| <b>Uniprot ID</b>        | <b>P47993</b>                            |
| <b>Sensitivity</b>       | 0.039 ng/ml                              |
| <b>Expiration Date</b>   | Please enquire.                          |