

Product Datasheet

Mouse pAMPK ELISA Kit (orb409238)

Description	Mouse Phosphorylated Adenosine Monophosphate Activated Protein Kinase, pAMPK ELISA Kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	pAMPK
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.