

Product Datasheet

Mouse PIK3CA ELISA Kit (orb409239)

Description	Mouse Phosphotylinosital 3 kinase,PI3K ELISA Kit
Reactivity	Mouse
Range	23.5 pg/mL-1500 pg/mL
Target	PIK3CA
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P42337
Sensitivity	5.8 pg/mL
Expiration Date	Please enquire.