

Product Datasheet

Mouse PCK1 ELISA Kit (orb409240)

Description	Mouse phosphoenolpyruvate carboxykinase,PCK ELISA Kit
Reactivity	Mouse
Range	1.56 pg/mL-100 pg/mL
Target	PCK1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9Z2V4
Sensitivity	0.39 pg/mL
Expiration Date	Please enquire.