

Product Datasheet

Mouse SEMA3A ELISA Kit (orb409267)

Description	Mouse Semaphorin-3A(SEMA3A) ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	SEMA3A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	O08665
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.