

## Product Datasheet

### Mouse Tissue factor ELISA Kit (orb409283)

|                          |                                                |
|--------------------------|------------------------------------------------|
| <b>Description</b>       | Mouse coagulation factor<br>III,FIII ELISA Kit |
| <b>Reactivity</b>        | Mouse                                          |
| <b>Range</b>             | 7.8 pg/mL-500 pg/mL                            |
| <b>Target</b>            | F3                                             |
| <b>Note</b>              | For research use only                          |
| <b>Application notes</b> | Detection Wavelength: 450<br>nm                |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates           |
| <b>Uniprot ID</b>        | <b>P20352</b>                                  |
| <b>Sensitivity</b>       | 1.95 pg/mL                                     |
| <b>Expiration Date</b>   | Please enquire.                                |