

Product Datasheet

Mouse GRP78 ELISA Kit (orb409300)

Description	Mouse Glucose Regulated Protein 78(GRP78)Elisa kit
Reactivity	Mouse
Range	0.45 ng/mL-30 ng/mL
Target	GRP78
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P20029
Sensitivity	0.11 ng/mL
Expiration Date	Please enquire.