

Product Datasheet

Mouse CXCL17 ELISA Kit (orb409314)

Description	Mouse chemokine (C-X-C motif) ligand 17 (CXCL17) ELISA kit
Reactivity	Mouse
Range	37.5 pg/mL-2400 pg/mL
Target	CXCL17
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q5UW37
Sensitivity	9.375 pg/mL
Expiration Date	Please enquire.