

Product Datasheet

Mouse HSP70 ELISA Kit (orb409322)

Description	Mouse Heat Shock Protein 70,Hsp-70 ELISA Kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	HSPA1A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q61696
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.