

Product Datasheet

Mouse HSP90AA1 ELISA Kit (orb409323)

Description	Mouse Heat Shock Protein 90,Hsp-90 ELISA Kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	HSP90AA1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P07901
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.