

Product Datasheet

Mouse Neurotrophin 3 ELISA Kit (orb409349)

Description	Mouse Neurotrophin 3,NT-3 ELISA kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	NT-3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P20181
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.