

Product Datasheet

Mouse HAVCR1 ELISA Kit (orb409355)

Description	Mouse Kidney injury molecule 1,Kim-1 ELISA Kit
Reactivity	Mouse
Range	0.9 ng/mL-60 ng/mL
Target	HAVCR1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, urine
Uniprot ID	Q5QNS5
Sensitivity	0.225 ng/mL
Expiration Date	Please enquire.