

Product Datasheet

Mouse GDF15 ELISA Kit (orb409356)

Description	Mouse Growth/differentiation factor 15(GDF15) ELISA kit
Reactivity	Mouse
Range	7.8 pg/ml -500 pg/ml
Target	GDF15
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9Z0J7
Sensitivity	1.95 pg/ml
Expiration Date	Please enquire.