

## Product Datasheet

### Mouse GH1 ELISA Kit (orb409358)

**Description** Mouse growth hormone,GH ELISA KIT

**Reactivity** Mouse

**Range** 3.12 pg/mL-200 pg/mL

**Target** GH1

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, cell culture  
supernates, cell lysates

**Uniprot ID** **P06880**

**Sensitivity** 0.78 pg/mL

**Expiration Date** Please enquire.