

## Product Datasheet

### Mouse PPIA ELISA Kit (orb409363)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | Mouse cyclophilin A,CyPA ELISA Kit                            |
| <b>Reactivity</b>        | Mouse                                                         |
| <b>Range</b>             | 15.6 pg/mL-1000 pg/mL                                         |
| <b>Target</b>            | PPIA                                                          |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                  |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates,<br>tissue homogenates |
| <b>Uniprot ID</b>        | <b>P17742</b>                                                 |
| <b>Sensitivity</b>       | 3.9 pg/mL                                                     |
| <b>Expiration Date</b>   | Please enquire.                                               |