

Product Datasheet

Mouse MYD88 ELISA Kit (orb409379)

Description	Mouse myeloid differentiation factor 88, MyD88 ELISA Kit
Reactivity	Mouse
Range	47 pg/mL-3000 pg/mL
Target	MYD88
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P22366
Sensitivity	11.7 pg/mL
Expiration Date	Please enquire.