

Product Datasheet

Mouse TREM2 ELISA Kit (orb409380)

Description	Mouse Triggering receptor expressed on myeloid cells 2(TREM2) ELISA kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	TREM2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q99NH8
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.