

## Product Datasheet

### Mouse AGER ELISA Kit (orb409412)

|                          |                                                                            |
|--------------------------|----------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Advanced glycosylation end product-specific receptor(AGER) ELISA kit |
| <b>Reactivity</b>        | Mouse                                                                      |
| <b>Range</b>             | 0.312 ng/ml-20 ng/ml                                                       |
| <b>Target</b>            | AGER                                                                       |
| <b>Note</b>              | For research use only                                                      |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                               |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                                          |
| <b>Uniprot ID</b>        | <b>Q62151</b>                                                              |
| <b>Sensitivity</b>       | 0.078 ng/ml                                                                |
| <b>Expiration Date</b>   | Please enquire.                                                            |