

Product Datasheet

Mouse CTLA4 ELISA Kit (orb409422)

Description	Mouse cytotoxic T lymphocyte associated antigen 4,CTLA-4 ELISA Kit
Reactivity	Mouse
Range	0.312 ng/mL-20 ng/mL
Target	CTLA4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P09793
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.