

Product Datasheet

Mouse ICAM1 ELISA Kit (orb409423)

Description	Mouse intercellular adhesion molecule 1,ICAM-1 ELISA KIT
Reactivity	Mouse
Range	78 ng/mL-5000 ng/mL
Target	ICAM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P13597
Sensitivity	19.5 ng/mL
Expiration Date	Please enquire.