

Product Datasheet

Mouse Cytochrome C ELISA Kit (orb409428)

Description	Mouse Cytochrome-C,Cyt-C ELISA Kit
Reactivity	Mouse
Range	34.3 pg/mL-2200 pg/mL
Target	CYCS
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P62897
Sensitivity	8.5 pg/mL
Expiration Date	Please enquire.