

Product Datasheet

Mouse SOCS-1 ELISA Kit (orb409430)

Description	Mouse suppressors of cytokine signaling 1(SOCS-1)ELISA Kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	SOCS1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O35716
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.