

Product Datasheet

Mouse FGL1 ELISA Kit (orb409436)

Description	Mouse Fibrinogen-like protein 1(FGL1) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	FGL1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q71KU9
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.