

## Product Datasheet

### Mouse SHBG ELISA Kit (orb409443)

|                          |                                                   |
|--------------------------|---------------------------------------------------|
| <b>Description</b>       | Mouse sex hormone-binding globulin,SHBG ELISA Kit |
| <b>Reactivity</b>        | Mouse                                             |
| <b>Range</b>             | 1.3 ng/mL-50 ng/mL                                |
| <b>Target</b>            | SHBG                                              |
| <b>Note</b>              | For research use only                             |
| <b>Application notes</b> | Detection Wavelength: 450 nm                      |
| <b>Sample Types</b>      | serum, plasma                                     |
| <b>Uniprot ID</b>        | <b>P97497</b>                                     |
| <b>Sensitivity</b>       | 0.33 ng/mL                                        |
| <b>Expiration Date</b>   | Please enquire.                                   |