

Product Datasheet

Mouse Vasoactive Intestinal Peptide ELISA Kit (orb409446)

Description	Mouse Vasoactive Intestinal Peptide,VIP ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Target	VIP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P32648
Sensitivity	1.95 pg/ml
Expiration Date	Please enquire.