

Product Datasheet

Mouse VEGFB ELISA Kit (orb409454)

Description	Mouse Vascular Endothelial cell Growth Factor B, VEGF-B ELISA kit
Reactivity	Mouse
Range	7.8 pg/mL-500 pg/mL
Target	VEGFB
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P49766
Sensitivity	0.195 pg/mL
Expiration Date	Please enquire.