

Product Datasheet

Mouse FIGF ELISA Kit (orb409456)

Description	Mouse Vascular Endothelial cell Growth Factor D,VEGF-D ELISA kit
Reactivity	Mouse
Range	0.312 pg/mL-20 pg/mL
Target	Vegfd
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P97946
Sensitivity	0.078 pg/mL
Expiration Date	Please enquire.