

Product Datasheet

Mouse KDR ELISA Kit (orb409458)

Description Mouse Vascular endothelial cell growth factor receptor 2, VEGFR-2/Flk-1 ELISA kit

Reactivity Mouse

Range 7.8 pg/mL-500 pg/mL

Target VEGFR-2/Flk-1

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma

Uniprot ID **P35918**

Sensitivity 1.95 pg/mL

Expiration Date Please enquire.