

Product Datasheet

Mouse FLT4 ELISA Kit (orb409459)

Description	Mouse Vascular endothelial cell growth factor receptor 3 (VEGFR-3/Flt-4) ELISA kit
Reactivity	Mouse
Range	78 pg/ml-5000 pg/ml.
Target	VEGFR-3/Flt-4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P35917
Sensitivity	19.5 pg/ml.
Expiration Date	Please enquire.