

Product Datasheet

Mouse ANGPTL4 ELISA Kit (orb409463)

Description	Mouse Angiopoietin-related protein 4(ANGPTL4) ELISA kit
Reactivity	Mouse
Range	3.12 pg/mL-200 pg/mL
Target	ANGPTL4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9Z1P8
Sensitivity	0.78 pg/mL
Expiration Date	Please enquire.