

## Product Datasheet

### Mouse THBS1 ELISA Kit (orb409475)

|                          |                                           |
|--------------------------|-------------------------------------------|
| <b>Description</b>       | Mouse thrombospondin 1,TSP-1<br>ELISA Kit |
| <b>Reactivity</b>        | Mouse                                     |
| <b>Range</b>             | 3.12 ng/ml-200 ng/ml                      |
| <b>Target</b>            | THBS1                                     |
| <b>Note</b>              | For research use only                     |
| <b>Application notes</b> | Detection Wavelength: 450 nm              |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates      |
| <b>Uniprot ID</b>        | <b>P35441</b>                             |
| <b>Sensitivity</b>       | 1.56 ng/ml                                |
| <b>Expiration Date</b>   | Please enquire.                           |