

Product Datasheet

Mouse PDGFC ELISA Kit (orb409481)

Description	Mouse platelet-derived growth factor CC (PDGF-CC) ELISA Kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	PDGFC
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q8CI19
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.