

Product Datasheet

Mouse PDGFRA ELISA Kit (orb409482)

Description	Mouse Platelet-Derived Growth Factor Soluble Receptor α , PDGFR- α ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	PDGFRA
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P26618
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.