

Product Datasheet

Mouse AMY2A ELISA Kit (orb409490)

Description	Mouse pancreatic alpha amylase (AMY2A)ELISA kit
Reactivity	Mouse
Range	9.38 mU/mL-600 mU/mL
Target	AMY2A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P00688
Sensitivity	2.3 mU/mL
Expiration Date	Please enquire.