

Product Datasheet

Mouse Insulin ELISA Kit (orb409492)

Description	Mouse Insulin,INS ELISA Kit
Reactivity	Mouse
Range	15.6 nIU/mL-1000 nIU/mL
Target	INS
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P01325
Sensitivity	3.9 nIU/mL
Expiration Date	Please enquire.