

Product Datasheet

Mouse IGF2 ELISA Kit (orb409494)

Description	Mouse Insulin-like growth factor 2,IGF-2 ELISA Kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	IGF2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P09535
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.