

Product Datasheet

Mouse IGFBP1 ELISA Kit (orb409495)

Description	Mouse Insulin-like growth factor binding protein 1,IGFBP-1 ELISA Kit
Reactivity	Mouse
Range	0.312 ng/mL-20 ng/mL
Target	IGFBP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P47876
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.