

Product Datasheet

Mouse IGFBP2 ELISA Kit (orb409496)

Description	Mouse Insulin-like growth factor binding protein 2,IGFBP-2 ELISA KIT
Reactivity	Mouse
Range	0.45 ng/mL-30 ng/mL
Target	IGFBP2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P47877
Sensitivity	0.11 ng/mL
Expiration Date	Please enquire.