

## Product Datasheet

### Mouse Pancreatic Lipase ELISA Kit (orb409506)

|                          |                                      |
|--------------------------|--------------------------------------|
| <b>Description</b>       | Mouse Pancreatic Lipase(PL)ELISA Kit |
| <b>Reactivity</b>        | Mouse                                |
| <b>Range</b>             | 62.5 mU/mL-4000 mU/mL                |
| <b>Target</b>            | PL                                   |
| <b>Note</b>              | For research use only                |
| <b>Application notes</b> | Detection Wavelength:<br>450 nm      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates    |
| <b>Uniprot ID</b>        | <b>Q6P8U6</b>                        |
| <b>Sensitivity</b>       | 15.6 mU/mL                           |
| <b>Expiration Date</b>   | Please enquire.                      |