

Product Datasheet

Mouse IDO ELISA Kit (orb409517)

Description	Mouse Indoleamine 2,3-dioxygenase 1(IDO1) ELISA kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	IDO1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P28776
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.