

Product Datasheet

Mouse iNOS ELISA Kit (orb409522)

Description	Mouse inducible nitric oxide synthase,iNOS ELISA KIT
Reactivity	Mouse
Range	7.8 IU/mL-500 IU/mL
Target	NOS2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P29477
Sensitivity	1.95 IU/mL
Expiration Date	Please enquire.