

Product Datasheet

Mouse PLA2G7 ELISA Kit (orb409537)

Description	Mouse lipoprotein-associated phospholipase A2,Lp-PLA2 ELISA Kit
Reactivity	Mouse
Range	7.8 pg/mL-500 pg/mL
Target	PLA2G7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q60963
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.