

Product Datasheet

Mouse NGAL ELISA Kit (orb409545)

Description	Mouse neutrophil gelatinase-associated lipocalin,NGAL ELISA Kit
Reactivity	Mouse
Range	12.5 pg/mL-800 pg/mL
Target	LCN2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, urine
Uniprot ID	P11672
Sensitivity	3.12 pg/mL
Expiration Date	Please enquire.