

Product Datasheet

Mouse TNF ELISA Kit (orb409547)

Description	Mouse Tumor necrosis factor α , TNF- α ELISA KIT
Reactivity	Mouse
Range	15.6 pg/ml-1000 pg/ml
Target	TNF
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P06804
Sensitivity	3.9 pg/ml
Expiration Date	Please enquire.