

Product Datasheet

Mouse TNFAIP6 ELISA Kit (orb409554)

Description	Mouse Tumor necrosis factor-inducible gene 6 protein(TNFAIP6) ELISA kit
Reactivity	Mouse
Range	62.5 pg/mL-4000 pg/mL
Target	TNFAIP6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	O08859
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.