

## Product Datasheet

### Mouse TGFA ELISA Kit (orb409556)

|                          |                                                                     |
|--------------------------|---------------------------------------------------------------------|
| <b>Description</b>       | Mouse transforming growth factor $\alpha$ , TGF- $\alpha$ ELISA Kit |
| <b>Reactivity</b>        | Mouse                                                               |
| <b>Range</b>             | 12.5 pg/mL-800 pg/mL                                                |
| <b>Target</b>            | TGFA                                                                |
| <b>Note</b>              | For research use only                                               |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                        |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates, tissue homogenates          |
| <b>Uniprot ID</b>        | <b>P48030</b>                                                       |
| <b>Sensitivity</b>       | 3.12 pg/mL                                                          |
| <b>Expiration Date</b>   | Please enquire.                                                     |