

Product Datasheet

Mouse TGF beta 1 ELISA Kit (orb409557)

Description	Mouse Transforming Growth factor β 1,TGF- β 1 ELISA kit
Reactivity	Mouse
Range	0.78 ng/mL-50 ng/mL
Target	TGFB1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates
Uniprot ID	P04202
Sensitivity	0.2 ng/mL
Expiration Date	Please enquire.